



Preparation of Nuclear Matrix on Slide for Immunofluorescence

- affix cells onto poly-L-lysine-coated coverslips as per standard protocol (e.g., cytospin)
- wash cells with PBS, 5 min
- incubate cells for 5 min with PBS containing 0.5% Triton X-100 (TX-100)
- remove the TX-100
- add ice-cold CSK buffer containing 0.5% TX-100 and incubate 5 min
- remove the CSK/T-100
- wash briefly twice with CSK buffer
- add DNase solution (1 mg/ml DNase I made in CSK buffer)
- incubate for 30 min
- remove DNase solution
- rinse and wash 2 x 5 min with PBS
- fix samples with 3% paraformaldehyde, 10 min
- continue with the regular immunofluorescence protocol

Solutions

- 1) CSK buffer: see protocols for buffers; add all protease inhibitors
- 2) DNase I solution: make up DNase I to 1 mg/ml in CSK buffer (will need ~250 μ l / coverslip)
- 3) all immunofluorescence solutions